



PROCAL
V5

PROCAL
DIRECT

RELEASE NOTES

VERSION 5.8.0.5

REVISION 0, SEPTEMBER 2019

RELEASE NOTES

Contents

Release Details:3

- Applicable Modules..... 3
- Modification Details..... 3
- Enhancements..... 3
- Corrections..... 4
- Obscure Bugs 4

Database4

- Enhancements..... 4
- Corrections..... 4

Documentation.....4

RELEASE NOTES

These notes describe in detail what elements were changed during the release process of applicable modules for the ProCalV5 Calibration and Maintenance and Quality Management Systems.

RELEASE DETAILS:

Applicable Modules

All application and system modules along with all documentation have been updated to version 5.8.0.4.

Modification Details

Enhancements

S-03525 Add the ability to create Instrument IDs from an Uploaded ad-hoc calibration on the Fluke Calibrator:

When Uploading an ad-hoc result and there is not a matching ID in the Company, you will have the option to create a new Instrument record along with the Calibration. The new Instrument will have the ID, Serial Number, (if entered on the Calibrator) and the Test Point Group.

S-06132 Show Equipment Description in dropdown lists on Item screens:

On the Instrument and Loop screens, the Equipment ID dropdown will now show the Equipment Description.

S-06133 Show Test Instrument Description on Calibration screen

On the Calibration screen, Test Instruments tab, the displayed Test Instrument list will now contain:

- Company
- Test Instrument ID
- Test Instrument Description

S-06134 Expand login screen's password field to more than 20 characters

To accommodate larger Passwords.

S-06135 Make list screens resizable and bigger

When on List screens you will now be able to expand the screen size.

S-06163 Add Shortcut keys for initiating a Calibration from a master record.

The shortcut keys of "CTRL+Q" will now start a Calibration for the master record that you are on.

S-06174 The display of all Logged on Users will now include the Time/Date of when they logged on.

This change will help admins when they need to Force a User off the system to free up a license.

RELEASE NOTES

Corrections

D-02295 If a Saved Query has the same name in multiple Companies and is Deleted in one Company it will be Deleted in all Companies

D-02287 Frequencies created in French that use 'Year' do not work

D-02302 "Send File via DeviceSync" file path fix

D-02336 Company screen, Options tab, the checkbox for "User Cannot Finalize Maint Results..." appears twice

D-02337 Cannot export via TMS grid if list contains a date in the selection formula.

S-06152 Export Audit Trail to Excel via the new grid (ProCal Direct)

S-06154 Task screen, opening the Task screen from the Work Order screen opens to the Notes tab

Obscure Bugs

There were no obscure bugs fixed in this release.

DATABASE

Enhancements

There were no enhancements made to the database in this release.

Corrections

There were no corrections made to the database for this release.

DOCUMENTATION

The User Manual was updated to include the additional features for this release